
pobm Documentation

Author

Oct 18, 2020

CONTENTS:

1	pobm package	1
1.1	Subpackages	1
1.1.1	pobm.obm package	1
1.1.1.1	Submodules	1
1.1.1.2	pobm.obm.burden	1
1.1.1.3	pobm.obm.complex	2
1.1.1.4	pobm.obm.desat	3
1.1.1.5	pobm.obm.general	4
1.1.1.6	pobm.obm.periodicity	5
1.1.1.7	Module contents	6
1.1.2	pobm.spo2 package	6
1.1.2.1	Submodules	6
1.1.2.2	pobm.spo2.single_biomarkers	6
1.1.2.3	Module contents	7
1.2	Submodules	7
1.3	pobm.main	7
1.4	pobm.prep	7
1.5	Module contents	8
2	Indices and tables	9
	Python Module Index	11
	Index	13

POBM PACKAGE

1.1 Subpackages

1.1.1 pobm.obm package

1.1.1.1 Submodules

1.1.1.2 pobm.obm.burden

class pobm.obm.burden.HypoxicBurdenMeasures (*begin*, *end*, *CT_Threshold=90*,
CA_Baseline=None)

Bases: object

Class that calculates Hypoxic Burden Features from spo2 time series. Suppose that the data has been preprocessed.

Parameters

- **begin** – List of indices of beginning of each desaturation event.
- **end** – List of indices of end of each desaturation event.
- **CT_Threshold** – Percentage of the time spent below the “CT_Threshold” % oxygen saturation level.
- **CA_Baseline** – Baseline to compute the CA feature. Default value is mean of the signal.

compute (*signal*)

Parameters **signal** – 1-d array, of shape (N,) where N is the length of the signal

Returns

HypoxicBurdenMeasuresResults class containing the following features:

- CA: Integral SpO2 below the xx SpO2 level normalized by the total recording time
- CT: Percentage of the time spent below the xx% oxygen saturation level
- POD: Percentage of oxygen desaturation events
- AODmax: The area under the oxygen desaturation event curve, using the maximum SpO2 value as baseline and normalized by the total recording time
- AOD100: Cumulative area of desaturations under the 100% SpO2 level as baseline and normalized by the total recording time

Example:

```
from pobm.obm.burden import HypoxicBurdenMeasures

# Initialize the class with the desired parameters
hypoxic_class = HypoxicBurdenMeasures(results_desat.begin, results_desat.end,
    ↪CT_Threshold=90, CA_Baseline=90)

# Compute the biomarkers
results_hypoxic = hypoxic_class.compute(spo2_signal)
```

1.1.1.3 pobm.obm.complex

```
class pobm.obm.complex.ComplexityMeasures (CTM_Threshold=0.25,      DFA_Window=20,
                                           M_Sampen=3,  R_Sampen=0.2,  M_ApEn=2,
                                           R_ApEn=0.25)
```

Bases: object

Class that calculates Complexity Features from spo2 time series. Suppose that the data has been preprocessed.

:param signal: 1-d array, of shape (N,) where N is the length of the signal CTM_Threshold: Radius of Central Tendency Measure. DFA_Window: Length of window to calculate DFA biomarker. M_Sampen: Embedding dimension to compute SampEn. R_Sampen: Tolerance to compute SampEn. M_ApEn: Embedding dimension to compute ApEn. R_ApEn: Tolerance to compute ApEn.

compute (signal) → pobm._ResultsClasses.ComplexityMeasuresResults

Parameters **signal** – 1-d array, of shape (N,) where N is the length of the signal

Returns

ComplexityMeasuresResults class containing the following features:

- ApEn: Approximate Entropy.
- LZ: Lempel-Ziv complexity.
- CTM: Central Tendency Measure.
- SampEn: Sample Entropy.
- DFA: Detrended Fluctuation Analysis.

Example:

```
from pobm.obm.complex import ComplexityMeasures

# Initialize the class with the desired parameters
complexity_class = ComplexityMeasures(CTM_Threshold=0.25, DFA_Window=20, M_
    ↪Sampen=3, R_Sampen=0.2, M_ApEn=2, R_ApEn=0.25)

# Compute the biomarkers
results_complexity = complexity_class.compute(spo2_signal)
```

1.1.1.4 pobm.obm.desat

class pobm.obm.desat.**DesaturationsMeasures** (*ODI_Threshold=3*)

Bases: object

Class that calculates the Desaturation Features from spo2 time series. Suppose that the data has been preprocessed.

Parameters **ODI_Threshold** – Threshold to compute Oxygen Desaturation Index.

compute (*signal*) → pobm._ResultsClasses.DesaturationsMeasuresResults

Parameters **signal** – 1-d array, of shape (N,) where N is the length of the signal

Returns

DesaturationsMeasuresResults class containing the following features:

- DL_u: Mean of desaturation length
- DL_sd: Standard deviation of desaturation length
- DA100_u: Mean of desaturation area using 100% as baseline.
- DA100_sd: Standard deviation of desaturation area using 100% as baseline
- DAmx_u: Mean of desaturation area using max value as baseline.
- DAmx_sd: Standard deviation of desaturation area using max value as baseline
- DD100_u: Mean of depth desaturation from 100%.
- DD100_sd: Standard deviation of depth desaturation from 100%.
- DDmx_u: Mean of depth desaturation from max value.
- DDmx_sd: Standard deviation of depth desaturation from max value.
- DS_u: Mean of the desaturation slope.
- DS_sd: Standard deviation of the desaturation slope.
- TD_u: Mean of time between two consecutive desaturation events.
- TD_sd: Standard deviation of time between 2 consecutive desaturation events.

Example:

```
from pobm.obm.desat import DesaturationsMeasures

# Initialize the class with the desired parameters
desat_class = DesaturationsMeasures(ODI_Threshold=3)

# Compute the biomarkers
results_desat = desat_class.compute(spo2_signal)
```

desaturation_detector (*signal*)

run desaturation detector, implemented by Dr. Joachim Behar

Parameters **signal** – The SpO2 signal, of shape (N,)

Returns

ODIMeasureResult class containing the following features:

- ODI: the average number of desaturation events per hour.
- begin: List of indices of beginning of each desaturation event.

- end: List of indices of end of each desaturation event.

`pobm.obm.desat.desat_embedding(begin, end)`

Help function for the class

Returns helper arrays containing the information about desaturation lengths and areas.

1.1.1.5 pobm.obm.general

class `pobm.obm.general.OverallGeneralMeasures` (*ZC_Baseline=None, percentile=1, M_Threshold=2, DI_Window=12*)

Bases: object

Class that calculates Overall General Features from spo2 time series. Suppose that the data has been preprocessed.

Parameters

- **ZC_Baseline** – Baseline for calculating number of zero-crossing points.
- **percentile** – Percentile to perform. For example, for percentile 1, the argument should be 1
- **M_Threshold** – Percentage of the signal M_Threshold % below median oxygen saturation. Typically use 1,2 or 5

compute (*signal*) → `pobm._ResultsClasses.OverallGeneralMeasuresResult`

Parameters **signal** – 1-d array, of shape (N,) where N is the length of the signal

Returns

OverallGeneralMeasuresResult class containing the following features:

- AV: Average of the signal.
- MED: Median of the signal.
- Min: Minimum value of the signal.
- SD: Std of the signal.
- RG: SpO2 range (difference between the max and min value).
- P: percentile.
- M: Percentage of the signal x% below median oxygen saturation.
- ZC: Number of zero-crossing points.
- DI: Delta Index.

Example:

```
from pobm.obm.general import OverallGeneralMeasures

# Initialize the class with the desired parameters
statistics_class = OverallGeneralMeasures(ZC_Baseline=90, percentile=1, M_
↳Threshold=2, DI_Window=12)

# Compute the biomarkers
results_statistics = statistics_class.compute(spo2_signal)
```


1.1.1.6 pobm.obm.periodicity

class pobm.obm.periodicity.PRSAMeasures (PRSA_Window=10, K_AC=2)

Bases: object

Function that calculates PRSA Features from spo2 time series. Suppose that the data has been preprocessed.

:param PRSA_Window: Fragment duration of PRSA. K_AC: Number of values to shift when computing autocorrelation

compute (signal) → pobm._ResultsClasses.PRSAResults

Parameters **signal** – 1-d array, of shape (N,) where N is the length of the signal

Returns

PRSAResults class containing the following features:

- PRSAc: PRSA capacity.
- PRSAad: PRSA amplitude difference.
- PRSAos: PRSA overall slope.
- PRSAsb: PRSA slope before the anchor point.
- PRSAsa: PRSA slope after the anchor point.
- AC: Autocorrelation.

Example:

```
from pobm.obm.periodicity import PRSAMeasures

# Initialize the class with the desired parameters
prsa_class = PRSAMeasures (PRSA_Window=10, K_AC=2)

# Compute the biomarkers
results_PRSA = prsa_class.compute(spo2_signal)
```

class pobm.obm.periodicity.PSDMeasures

Bases: object

Function that calculates PSD Features from spo2 time series. Suppose that the data has been preprocessed.

compute (signal) → pobm._ResultsClasses.PSDResults

:param signal: The SpO2 signal, of shape (N,)

Returns

PSDResults class containing the following features:

- PSD_total: The amplitude of the spectral signal.
- PSD_band: The amplitude of the signal multiplied by a band-pass filter between 0.014 and 0.033 Hz.
- PSD_ratio: The ratio between PSD_total and PSD_band.
- PDS_peak: The max value of the FFT into the band 0.014-0.033 Hz.

Example:

```
from pobm.obm.periodicity import PSDMeasures

# Initialize the class with the desired parameters
psd_class = PSDMeasures()

# Compute the biomarkers
results_PSD = psd_class.compute(spo2_signal)
```

1.1.1.7 Module contents

1.1.2 pobm.spo2 package

1.1.2.1 Submodules

1.1.2.2 pobm.spo2.single_biomarkers

`pobm.spo2.single_biomarkers.apen` (*signal*, *M_ApEn*=2, *R_ApEn*=0.25)

Compute the approximate entropy, according to the paper Utility of Approximate Entropy From Overnight Pulse Oximetry Data in the Diagnosis of the Obstructive Sleep Apnea Syndrome

Parameters *signal* – 1-d array, of shape (N,) where N is the length of the signal

Returns ApEn

`pobm.spo2.single_biomarkers.dfa` (*signal*, *DFA_Window*=20)

Compute DFA

Parameters

- **signal** – 1-d array, of shape (N,) where N is the length of the signal
- **DFA_Window** – Length of window to calculate DFA biomarker.

Returns DFA

`pobm.spo2.single_biomarkers.lempel_ziv` (*signal*)

Compute lempel-ziv, according to the paper Non-linear characteristics of blood oxygen saturation from nocturnal oximetry for obstructive sleep apnoea detection

Parameters *signal* – 1-d array, of shape (N,) where N is the length of the signal

Returns LZ

`pobm.spo2.single_biomarkers.odi` (*signal*, *ODI_Threshold*=3)

Calculates the ODI from spo2 time series. Suppose that the data has been preprocessed.

Parameters

- **signal** – The SpO2 signal, of shape (N,)
- **ODI_Threshold** – Threshold to compute Oxygen Desaturation Index.

Returns ODI

`pobm.spo2.single_biomarkers.sampen` (*signal*, *M_Sampen*=3, *R_Sampen*=0.2)

Compute the Sample Entropy

Parameters

- **signal** – 1-d array, of shape (N,) where N is the length of the signal
- **M_Sampen** – Embedding dimension to compute SampEn.

- **R_Sampen** – Tolerance to compute SampEn.

Returns SampEn

1.1.2.3 Module contents

1.2 Submodules

1.3 pobm.main

1.4 pobm.prep

`pobm.prep.block_data` (*signal*, *threshold=50*)

Apply a block data filter to the SpO2 signal.

Parameters

- **signal** – 1-d array, of shape (N,) where N is the length of the signal
- **(Optional)** (*threshold*) – threshold parameter for block data filter.

Returns preprocessed signal, 1-d numpy array.

`pobm.prep.dfilter` (*signal*, *Diff=4*)

Apply Delta Filter to the signal.

Parameters

- **signal** – 1-d array, of shape (N,) where N is the length of the signal
- **Diff** – parameter of the delta filter.

Returns preprocessed signal, 1-d numpy array.

`pobm.prep.median_spo2` (*signal_spo2*, *FilterLength=9*)

Apply a median filter to the SpO2 signal. Median filter used to smooth the spo2 time series and avoid sporadic increase/decrease of spo2 which could affect the detection of the desaturations. Assumption: any missing/abnormal values are represented as 'np.nan'

Parameters

- **signal** – 1-d array, of shape (N,) where N is the length of the signal
- **(Optional)** (*FilterLength*) – The length of the filter.

Returns preprocessed signal, 1-d numpy array.

`pobm.prep.resamp_spo2` (*signal*, *OriginalFreq*)

Resample the SpO2 signal to 1Hz. Assumption: any missing/abnormal values are represented as 'np.nan'

Parameters

- **signal** – 1-d array, of shape (N,) where N is the length of the signal
- **OriginalFreq** – the original frequency.

Returns resampled signal, 1-d numpy array, the resampled spo2 time series at 1Hz

`pobm.prep.set_range` (*signal*, *Range_min=50*, *Range_max=100*)

Range function. Remove values lower than 50 or greater than 100, considered as non-physiological

Parameters **signal** – 1-d array, of shape (N,) where N is the length of the signal

Returns preprocessed signal, 1-d numpy array.

1.5 Module contents

INDICES AND TABLES

- `genindex`
- `modindex`
- `search`

PYTHON MODULE INDEX

p

- pobm, 8
- pobm.obm, 6
- pobm.obm.burden, 1
- pobm.obm.complex, 2
- pobm.obm.desat, 3
- pobm.obm.general, 4
- pobm.obm.periodicity, 5
- pobm.prep, 7
- pobm.spo2, 7
- pobm.spo2.single_biomarkers, 6

A

`apen()` (in module *pobm.spo2.single_biomarkers*), 6

B

`block_data()` (in module *pobm.prep*), 7

C

`ComplexityMeasures` (class in *pobm.obm.complex*), 2

`compute()` (*pobm.obm.burden.HypoxicBurdenMeasures* method), 1

`compute()` (*pobm.obm.complex.ComplexityMeasures* method), 2

`compute()` (*pobm.obm.desat.DesaturationsMeasures* method), 3

`compute()` (*pobm.obm.general.OverallGeneralMeasures* method), 4

`compute()` (*pobm.obm.periodicity.PRSAMeasures* method), 5

`compute()` (*pobm.obm.periodicity.PSDMeasures* method), 5

D

`desat_embedding()` (in module *pobm.obm.desat*), 4

`desaturation_detector()` (*pobm.obm.desat.DesaturationsMeasures* method), 3

`DesaturationsMeasures` (class in *pobm.obm.desat*), 3

`dfa()` (in module *pobm.spo2.single_biomarkers*), 6

`dfilter()` (in module *pobm.prep*), 7

H

`HypoxicBurdenMeasures` (class in *pobm.obm.burden*), 1

L

`lempel_ziv()` (in module *pobm.spo2.single_biomarkers*), 6

M

`median_spo2()` (in module *pobm.prep*), 7

module

pobm, 8
pobm.obm, 6
pobm.obm.burden, 1
pobm.obm.complex, 2
pobm.obm.desat, 3
pobm.obm.general, 4
pobm.obm.periodicity, 5
pobm.prep, 7
pobm.spo2, 7
pobm.spo2.single_biomarkers, 6

O

`odi()` (in module *pobm.spo2.single_biomarkers*), 6
`OverallGeneralMeasures` (class in *pobm.obm.general*), 4

P

pobm
 module, 8
pobm.obm
 module, 6
pobm.obm.burden
 module, 1
pobm.obm.complex
 module, 2
pobm.obm.desat
 module, 3
pobm.obm.general
 module, 4
pobm.obm.periodicity
 module, 5
pobm.prep
 module, 7
pobm.spo2
 module, 7
pobm.spo2.single_biomarkers
 module, 6
`PRSAMeasures` (class in *pobm.obm.periodicity*), 5
`PSDMeasures` (class in *pobm.obm.periodicity*), 5

R

`resamp_spo2()` (*in module pobm.prep*), [7](#)

S

`sampen()` (*in module pobm.spo2.single_biomarkers*), [6](#)

`set_range()` (*in module pobm.prep*), [7](#)